Kimtech™ G3 Sterile White Nitril- Handschuhe

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Chemical Permeation Table (printed 2021-02-15)
Standard [EN16523-1]

Glove	Concentration	CAS Number	Kimtech™ G3 Sterile White Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Acetic Acid	99%	64-19-7	5,0	94,0%
Ammonium Hydroxide	25%	1336-21-6	16	2,7%
Chloroform	99%	67-66-3	<1	92,0%
Citric Acid	30%	77-92-9	>480	-18,0%
Cyclohexane	99%	100-82-7	55	52,0%
Dimethyl Sulphoxide	99,50%	67-68-5	11	69,0%
Ethanol	70%	64-17-5	32	56,0%

Glove	Concentration	CAS Number	Kimtech™ G3 Sterile White Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Ethidium Bromide	1%	1239-45-8	>480	15,0%
Ethyl Acetate	99%	141-78-6	1,1	76,0%
Formaldehyde	37%	50-00-0	>480	29,0%
Hydrochloric Acid	30%	7647-01-0	>480	22,0%
Hydrogen Peroxide	30%	7722-84-1	>480	34,0%
Isopropanol	70%	67-63-0	110	49,0%
Isopropanol	99%	67-63-0	34	59,0%
Methanol	99%	67-56-1	4,6	75,0%
Nitric Acid	65%	7697-37-2	7,1	95,0%
Perchloric Acid	70%	7601-90-3	>480	9,9%
Sodium Hydroxide	40%	1310-73-2	>480	-13,0%

Glove	Concentration	CAS Number	Kimtech™ G3 Sterile White Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Sodium Hypochlorite	14%	7681-52-9	>480	13,0%
Sulphuric Acid	50%	7664-93-9	>480	5,2%
Toluene	99%	108-88-3	<1	77,0%
n-Heptane	99%	142-82-5	43	41,0%

Disclaimer: All data provided is based on results of tests performed in accordance with the relevant test standard (Chemical Permeation: EN16523-1 or EN374-3; Degradation: EN374-4), by an independent laboratory which has approval from a notified body under the CE Regulation (Or Directive) for Personal Protective Equipment. These tests may not adequately replicate any specific conditions of use, and because KIMBERLY-CLARK PROFESSIONAL™ has no detailed knowledge or control over the conditions of end use, any of the data provided must be considered on an advisory basis only, and KIMBERLY-CLARK PROFESSIONAL* must decline any liability.

Ref: Cytotoxic Drug Permeation testing of Kimtech™ gloves


Dear Valued Customer,

Thank you for your enquiry about our gloves and their use for protection against chemical splash hazards by cytotoxic drugs.

We must caution that the selection of the most appropriate glove for a specific task should be carried out by a trained safety professional following a full risk assessment. It is Kimberly-Clark's intent to provide data which can enable trained professionals to make an informed choice.

ASTM D 6978-05 (Standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs) was utilised to complete the chemotherapy drug performance assessment on the above glove products. Testing was performed on the cuff area of the gloves, under conditions of continuous contact. UV/VIS Spectrometry was used to measure the absorbance of the challenge chemicals through the specimens, into the collection medium.

K-C Confidential. The information contained in this document is intended for the use of the addressee(s) only and may contain privileged, confidential, or proprietary information that should not be distributed publically without Kimberly-Clark's written consent.

Results are summarised below:

Kimtech™G3 Sterile White Nitrile gloves

Chemotherapy Drug	Concentration (mg/ml)	Breakthrough Time (mins)
Carmustine	3.3	30.4
Cisplatin	1.0	No breakthrough up to 240 mins
Cyclophosphamide (Cytoxan)	20.0	No breakthrough up to 240 mins
Dacarbazine (DTIC)	10.0	No breakthrough up to 240 mins
Doxorubicin Hydrochloride	2.0	No breakthrough up to 240 mins
Etoposide (Toposar)	20.0	No breakthrough up to 240 mins
Fluorouracil	50.0	No breakthrough up to 240 mins
Ifosfamide	50.0	No breakthrough up to 240 mins
Mitoxantrone	2.0	No breakthrough up to 240 mins
Paclitaxel (Taxol)	6.0	No breakthrough up to 240 mins
Thiotepa	10.0	No breakthrough up to 240 mins
Vincristine Sulphate	1.0	No breakthrough up to 240 mins

Tests were carried out by an independent laboratory, under laboratory test conditions. These tests may not adequately replicate any specific condition of use. As KIMBERLY-CLARK* has no detailed knowledge or control over the conditions of end use, this data must be considered advisory only and KIMBERLY-CLARK* must decline any liability.

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¹ Nitril ist ein synthetisches Material, das viele Eigenschaften mit Naturkautschuk-Latex gemein hat, sich aber von diesem durch mehrere signifikante Vorteile unterscheidet: hoher Tragekomfort, gute Stichfestigkeit, hohe Abriebfestigkeit ohne Beeinträchtigung der Tastempfindlichkeit oder der elektrostatisch dissipativen Eigenschaften. ² AQL-Wert gemäß Definition nach ISO 2859-1 für Attributprüfung anhand von Proben. ©/™ Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. © KCWW. Publication code: ID4410.04 DE 11.19
Kimberly-Clark Europe Ltd.
Kimberly-Clark Europe Limited
TÜV SÜD Product Service GmbH Zertifizierstellen (0123)